

01/2008:20403
corrected 8.0

2.4.3. CALCIUM

All solutions used for this test are prepared with distilled water R.

To 0.2 mL of *alcoholic calcium standard solution* (100 ppm Ca) R add 1 mL of *ammonium oxalate solution* R. After 1 min add a mixture of 1 mL of *dilute acetic acid* R and 15 mL of the prescribed solution or of a solution containing the prescribed quantity of the substance to be examined, and shake. Prepare a standard in the same manner using a mixture of 10 mL of *aqueous calcium standard solution* (10 ppm Ca) R, 1 mL of *dilute acetic acid* R and 5 mL of *distilled water* R.

After 15 min, any opalescence in the test solution is not more intense than that in the standard.

01/2008:20404



2.4.4. CHLORIDES

To 15 mL of the prescribed solution add 1 mL of *dilute nitric acid* R and pour the mixture as a single addition into a test-tube containing 1 mL of *silver nitrate solution* R2. Prepare a standard in the same manner using 10 mL of *chloride standard solution* (5 ppm Cl) R and 5 mL of *water* R. Examine the tubes laterally against a black background.

After standing for 5 min protected from light, any opalescence in the test solution is not more intense than that in the standard.

01/2008:20405



2.4.5. FLUORIDES

Introduce into the inner tube of the apparatus (see Figure 2.4.5.-1) the prescribed quantity of the substance to be examined, 0.1 g of acid-washed *sand* R and 20 mL of a mixture of equal volumes of *sulfuric acid* R and *water* R. Heat the jacket containing *tetrachloroethane* R maintained at its boiling point (146 °C). Heat the steam generator and distil, collecting the distillate in a 100 mL volumetric flask containing 0.3 mL of 0.1 M *sodium hydroxide* and 0.1 mL of *phenolphthalein solution* R. Maintain a constant volume (20 mL) in the tube during distillation and ensure that the distillate remains alkaline, adding 0.1 M *sodium hydroxide* if necessary. Dilute the distillate to 100 mL with *water* R (test solution). Prepare a standard in the same manner by distillation, using 5 mL of *fluoride standard solution* (10 ppm F) R instead of the substance to be examined. Into two glass-stoppered cylinders introduce 20 mL of the test solution and 20 mL of the standard and 5 mL of *aminomethylazarindiacetic acid reagent* R.

After 20 min, any blue colour in the test solution (originally red) is not more intense than that in the standard.




Figure 2.4.5.-1. – Apparatus for limit test for fluorides
Dimensions in millimetres

01/2008:20406



2.4.6. MAGNESIUM

To 10 mL of the prescribed solution add 0.1 g of *disodium tetraborate* R. Adjust the solution, if necessary, to pH 8.8 to pH 9.2 using *dilute hydrochloric acid* R or *dilute sodium hydroxide solution* R. Shake with 2 quantities, each of 5 mL, of a 1 g/L solution of *hydroxyquinoline* R in *chloroform* R, for 1 min each time. Allow to stand. Separate and discard the organic layer. To the aqueous solution add 0.4 mL of