

1. 약사법 제31조제1항의 규정에 의한 제조·품질관리에 필요한 시설에 관한 자료

-Site master file을 첨부하였습니다.

- 제조사인 Jing Jing Pharmaceutical Co., Ltd에서 Restricted Part에 제조소 평면도, 작업환경 관리구역 표시도면, 공조시설계통도, 압축공기계통도 및 용수처리계통도를 식약처로 직접 송부할 예정입니다.

SungWun Pharmacopia

2. 품목별로 실시상황이 별표 1의2 의약품제조 및 품질관리기준에 적합하거나 이와 동등이상임을

입증하는 자료

해당 자료는 민원 신청 후 식약처 홈페이지를 통해 제조사가 직접 업로드 할 예정입니다.

SungWun Pharmacopia

3. 물리화학적 특성과 안정성자료

3.2.S.1 일반정보

엘-오르니틴-엘-아스파르트산은 간 보호 및 해독작용을 하여 간경변 및 간성뇌병증 치료, 만성감염 해독의 보조 치료제 등으로 사용되며, Jing Jing Pharmaceutical Co., Ltd에서 연구 개발한 제품이다.

3.2.S.1.1 명칭

제품명 : 엘-오르니틴-엘-아스파르트산


기타 이름 : L-ornithine-L-aspartate

화학명 : (2S)-2-aminobutanedioic acid;(2S)-2,5-diaminopentanoic acid

CAS No. : 3230-94-2

3.2.S.1.2 구조

- 구조식



▪ Molecular Formula: $C_9H_{19}N_3O_6$

▪ Molecular Weight :265.26

3.3.S.1.3 일반적 특성

엘-오르니틴-엘-아스파르트산의 구조는 원소분석 데이터를 기본으로 하여 핵자기공명스펙트럼 (NMR), 질량분석스펙트럼(MS), 원소분석, X선 회절분석(XRD), 적외선스펙트럼(IR)을 종합하여 화합물의 구조를 규명하였다.

SungWun Pharmacopia

3.2.S.3 특성

3.2.S.3.1 구조 및 기타 특성

1) 핵자기공명스펙트럼 (NMR)

엘-오르니틴-엘-아스파르트산을 AV II- 500 BKUKER nuclear magnetic resonance spectrometer를 사용하여 ^1H NMR, ^{13}C NMR을 측정하였다.

- ^1H NMR data

Chemical shift(δ)		Proton number	Multiplicity	Hydrogen identification
Sample (C552202010))	Standard product (07125-SLBW4848)			
1.716-1.743	1.694-1.709	2	m	H-2
1.787-1.805	1.769-1.783	2	m	H-3
2.622-2.809	2.615-2.783	2	m	H-8
3.011-3.041	2.989-3.019	2	t	H-1
3.745-3.769	3.727-3.751	1	t	H-4
3.858-3.883	3.838-3.862	1	dd	H-7

- ^{13}C NMR data

Chemical shift(δ)		Carbon type	Carbon identification
Sample (C552202010))	Standard product (07125-SLBW4848)		
22.12	22.13	Secondary carbon	C-2
26.79	26.81	Secondary carbon	C-3
35.93	35.93	Secondary carbon	C-1
38.24	38.27	Secondary carbon	C-8
51.56	51.62	Tertiary carbon	C-7
53.48	53.51	Tertiary carbon	C-4
173.53	173.51	Quaternary carbon	C-5
173.53	173.64	Quaternary carbon	C-6
176.92	176.96	Quaternary carbon	C-9

2) 질량분석스펙트럼 (Mass spectrum)

Thermofisher LTQ를 사용하여 질량 분석하였다.

Batch number	Mass-to-charge ratio(m/z)	Relative abundance	Remarks
Sample	264.90	100	---
Standard product	264.18	100	---

3) 원소분석

엘-오르니틴-엘-아스파르트산의 원소분석을 실시하고, 이론값과 측정값을 다음의 표로 정리하였다.

그 결과 시료와 대조품 C, H, N의 비율이 기준과 일치함을 확인하였다.

Element		Actual test value (%)	Theoretical value (%)	Atomic weight
Standard	C	39.49%	40.71%	12.0107
	H	7.50%	7.22%	1.00794
	N	15.21%	15.84%	15.9999
sample	C	39.96%	40.71%	12.0107
	H	7.49%	7.22%	1.00794
	N	15.63%	15.84%	15.9999

4) X선 회절분석(X-ray Diffraction)




분석 결과, 시료의 X-선 회절 피크의 위치가 기준값과 동일하며, 결정의 경우도 비슷함을 보였다.

5) 적외부스펙트럼(IR)

엘-오르니틴-엘-아스파르트산의 적외선 실험 결과, 시료가 기준물질의 적외선 스펙트럼과 일치함을 확인하였다.

7) 유연물질


엘-오르니틴-엘-아스파르트산에 대해 시험한 유연물질 결과를 다음과 같이 정리하였다.

C552202010		C552202011		C552202012	
					
sample solution	standard solution	sample solution	standard solution	sample solution	standard solution

- 결과

Maximum Daily Dose	Reporting Threshold	Identification Threshold	Qualification Threshold
≤ 2g/day	0.05%	0.1% or 1.0mg Per day intake (whichever is lower)	0.15% or 1.0mg Per day intake (whichever is lower)
> 2g/day	0.03%	0.05%	0.05%

[ICH Q3AR2, Attachment 1 Thresholds]




品名: 11-庚基-5-氧-2-萘酚
批号: (337)02010
检验人: 闫雅菲
检验日期: 2022.02.11

2.5% 0.05% 0.025%

2ml 1ml

sample solution Standard solution (1ml)



品名: 11-庚基-5-氧-2-萘酚
批号: (337)02010
检验人: 闫雅菲
检验日期: 2022.02.11

2.5% 0.025% 0.015%

1ml 0.6ml

sample solution standard solution (1ml) standard solution (0.6ml)

4. 제조방법, 포장, 용기, 취급상의 주의사항 등에 관한 자료


3.2.S.2. 제조

3.2.S.2.1 제조원

Jing Jing Pharmaceutical Co., Ltd


3.2.S.2.2 제조공정 및 공정관리

1) 제조공정 기술



2) 합성 경로

엘-오르니틴-엘-아스파르트산의 합성 과정은 다음과 같다.



3.2.S.2.3 원료관리

Restricted Part로 제조사 (Jing Jing Pharmaceutical Co., Ltd)에서 식약처 전자메일로 송부할 예정입니다.

3.2.S.2.4 주요공정 및 중간체 관리

Restricted Part로 제조사(Jing Jing Pharmaceutical Co., Ltd)에서 식약처 전자메일로 송부할 예정입니다.

3.2.S.2.5. 공정 밸리데이션 및 평가

Restricted Part로 제조사(Jing Jing Pharmaceutical Co., Ltd)에서 식약처 전자메일로 송부할 예정입니다.

3.2.S.2.6 제조공정 개발

제조공정은 PV 배치 생산을 통해 성공적으로 재현했으며 독일약전(DAB) 기준에 적합한 제품을 생산할 수 있다.

3.2.S.6. 용기 및 포장

1) 포장정보

엘-오르니틴-엘-아스파르트산은 polyethylene 백에 1차 포장 후 fiber-drum에 외포장하여 출하된다.

- Specification of Inner Packing Material

Item	Standard
Appearance	The surface should be smooth and uniform in color, and there should be no perforation, foreign matter, odor, or adhesion. The heat-sealing part of the bag should be flat and free of false sealing.
identify	should be consistent with the control pattern
Specifications Dimensions	Lenght(mm) 900 Width(mm) 600 Thickness (mm) 0.08
Microbial limit	The total number of aerobic bacteria ≤ 1000 cfu/ml The total number of mold and yeast ≤ 100 cfu/ml

- Specification of Fiber drum

Test item	Specification
Appearance	The paper barrel should be round and free of defects and cracks such as obvious out-of-roundness, concave deflation, skew, etc. The barrel body is smooth, no mechanical damage, no wrinkle, no glue opening. Paint spreads evenly. No paint leakage, no bubbles, no obvious sagging. Round curled edges without paper tongue. Closing and rear lid and barrel body are well sealed, and the barrel body is rolled with imported cardboard paper. The inside and outside of the barrel should be clean, free of impurities and oil stains.
Dimensions	Barrel body: 400mm \pm 3mm Inner height: 550mm \pm 3mm Outer height: 570mm \pm 4mm Weight: 3.15 \pm 0.3kg

5. 원료의약품의 시험성적서, 분석방법, 사용된 용매 등에 관한 자료

3.2.S.4. 원료의약품의 관리

3.2.S.4.1 기준

엘-오르니틴-엘-아스파르트산은 DAB 기준을 따른다.

Test		Specification	Methods
Appearance		White crystal or crystalline powder	DAB
Identification	A. Optical Rotation, °	+26.5 ~ +29.0 °	Ph. Eur. 2.2.7
	B. IR	corresponds to standard	Ph. Eur. 2.2.24
	C. Ninhydrin reaction	Violet	DAB
	D. Mercuric acetate reaction	White Precipitation	DAB
	E. Molybdophosphoric acid reaction	Yellow Precipitation	DAB
Purity	1) Clarity / Coloration	Clear, Colourless	Ph. Eur. 2.2.1, 2.2.2(II)
	2) pH	6.0 ~ 7.0	Ph. Eur. 2.2.3
	3) Related substance (TLC)	Spots other than the main spots in the test solution are not larger than the spots obtained in the standard solution	Ph. Eur. 2.2.27
	4) Chloride	≤ 300 ppm	Ph. Eur. 2.2.4
	5) Sulfate	≤ 200 ppm	Ph. Eur. 2.4.13
	6) Ammonium	≤ 400 ppm	DAB N 2.4.1
	7) Iron	≤ 30 ppm	Ph. Eur. 2.4.9
	8) Heavy metals	≤ 10 ppm	Ph. Eur. 2.4.8
	9) Water	≤ 7.0 %	Ph. Eur. 2.5.12
	10) Sulfated ash	≤ 0.2 %	Ph. Eur. 2.4.14
Assay		98.0 ~ 102.0 %	Ph. Eur. 2.2.20
Residual solvent(Methanol)		≤ 3,000 ppm	In-house

3.2.S.4.2. 시험방법

Gas chromatography		Condition
Detector		Flame ionization(FID)
Column		DB-624 Capillary(G43), 0.25 mm x 30 m, 1.4 μ m
Temp.	Detector	280 $^{\circ}$ C
	Oven	40 $^{\circ}$ C(5 min) \rightarrow 240 $^{\circ}$ C(Rate 20 $^{\circ}$ C/ min)(10 min)
Carrier gas		nitrogen
Split Ratio		20 : 1
Flows		1) Hydrogen : 40 mL/min, 2) Air : 400.0 mL/min

Head-space	Condition
Sample equilibration temperature	80 $^{\circ}$ C
Quantitative loop temperature	120 $^{\circ}$ C
Transmission line temperature	120 $^{\circ}$ C
Equilibrating Time	20 min
Injection Time	1 min(mL)
GC cycle time	25 min

•Calculation

$$\text{Residue of MeOH (ppm)} = \frac{A_t \times W_s \times 1 \times 10^6}{A_s \times 1000 \times W_t} \times \frac{P_s}{100}$$

A_t : Peak area of residual solvent in the sample solution

A_s : Peak area of residual solvent in standard solution

1 : Dilution factor of sample

1000 : Dilution of standard solution

W_t : Weight of Sample (mg)

W_s : Weight of Standard (mg)

P_s : Purity of Standard

3.2.S.4.3 시험방법 밸리데이션

분석방법은 DAB에 준거한다.

Test Item	Criteria	Validation results	Result
System suitability test	The peak area measurement value RSD% (n=6) of repeated injection shall not exceed 10.0%; the number of theoretical plates (N) \geq 3000	Peak area repeatability RSD of 1.1% The minimum number of theoretical plates is 17425	Conform
Specificity	The blank solvent has no interference with methanol detection, and no other impurity peaks in the test solution interfere with each known impurity peak.	no interference	Conform
Detection of Limit	S/N is about 3, which can meet the testing requirements	The detection limit is 0.6010 μ g/ml, which is equivalent to 3ppm of the test sample	Conform
Detection of Quantitation	S/N is about 10, which can meet the testing requirements; Take the limit of quantification solution for 6 consecutive injections, and calculate the RSD of the peak area \leq 10.0% and the RSD of the retention time \leq 10.0%.	The limit of quantification is 1.8030 μ g/ml, which is equivalent to the percentage content of the test product of 9ppm The RSD for the limit of quantification precision was 1.8% RSD for retention time is 0.1%	Conform
Linearity and range	$R \leq 0.990$; The Y-axis intercept is within 25% of the 100% response value; Response factor RSD \leq 10%	Linear equation: $y = 846.56x - 7.1573$ R is 0.99996 The percentage of Y-intercept to 100% response value is 0.870% Linear range 0.030080~1.20320mg/ml Response factor RSD of 4.4.	Conform
Precision	Repeatability: The RSD of the peak area for 6 consecutive injections of the reference solution is \leq 10.0%; the RSD of the impurity content of the 6 samples of the test solution is \leq 10.0%. Intermediate precision: Impurity content of the test solution	Repeatability: The RSD of the reference solution is 1.1% The RSD of the test solution is 0.7% Intermediate precision: The RSD of the test solution (n=6) is 1.2% The RSD of the test solution (n=12) is	Conform

	RSD% (n=6) \leq 10.0%, RSD% (n=12) \leq 15.0%	1.0%	
Accuracy test	According to the external standard method, the detected amount and recovery rate of each impurity were calculated. The average recovery rate is between 80% and 120%, and the relative standard deviation should not exceed 10.0%	50% recovery was 102%, RSD (n=3) of 1.0% 100% recovery was 99%, RSD (n=3) of 0.6% 150% recovery was 99%, RSD (n=3) of 1.6% The average recovery was 100%, and the RSD (n=9) was 2.0%	Conform
Durability	Under each condition, the RD of methanol content in the solution of the spiked test sample shall not exceed 15.0%)	When there is a slight change in the measurement conditions, the theoretical plate number (N) of the reference substance is more than 3000, and the RD is less than 10.0%, all of which meet the requirements, and the measurement results are within the acceptable range. Including: different flow rates (1.315~1.715ml/min), different detector temperatures (275~285 $^{\circ}$ C), different inlet temperatures (115~125 $^{\circ}$ C), different headspace equilibration times (18~226min), Different headspace temperature (78 ~ 82 $^{\circ}$ C).	Conform

3.2.S.4.4 배치 분석

1) Lot COA:

	Batch Number	Batch Size	Manufacture Date	Analysis Date
PV1	C552202010	400.84 kg	2022.02.08	2022.02.10~2022.02.20
PV2	C552202011	400.84 kg	2022.02.09	2022.02.10~2022.02.20
PV3	C552202012	400.84 kg	2022.02.09	2022.02.11~2022.02.20

2) Result of Test

Batch number		C552202010	C552202011	C552202012
Manufacturing date		2022.02.08	2022.02.09	2022.02.09
Batch quantity		400.84kg	400.84kg	400.84kg
Testing Item	Standad	Result	Result	Result
Description	White crystal or crystalline powder	White crystalline powder	White crystalline powder	White crystalline powder
pH	6.0 ~ 7.0	6.3	6.3	6.3
Identifi- cation	IR	should be consistent with the standard infrared spectrum	Consistent with standard infrared spectrum	Consistent with standard infrared spectrum
	Optical Rotation, $^{\circ}$	+26.5 ~ +29.0	28.0	28.3
Clarity of the solution	≤No. 1 Turbidity Standard	<No. 1 Turbidity Standard	<No. 1 Turbidity Standard	<No. 1 Turbidity Standard
color of solution	≤ B9	< B9	< B9	< B9
Wate %	≤ 7.0	0.81	0.66	0.77

Residue on ignition	≤ 0.2	0.02	0.03	0.01
Chloride,ppm	≤ 300	< 300	< 300	< 300
Sulfate(ppm)	≤ 200	< 200	< 200	< 200
Ammonium ,ppm	≤ 400	< 400	< 400	< 400
Heavy metal,ppm	≤ 10	< 10	< 10	< 10
Iron, ppm	≤ 30	< 30	< 30	< 30
Other Amino Acids,%	Spots other than the main spots in the test solution are not larger than the spots obtained in the standard solution	Compliance	Compliance	Compliance
Residual solvent methanol ppm	≤3000	91	95	89
Assay, %	98.0 ~ 102.0	98.6	98.9	99.0

3.2.S.4.5 기준 설정근거

엘-오르니틴-엘-아스파르트산은 DAB 기준을 따르며, 잔류용매는 ICH Q3C를 따라 진행하였다.
일반 시험법은 DAB 및 EP을 따라 진행하였다.

SungWun Pharmacopia

3.2.S.5 표준품 또는 표준물질

Standard substance	Source	Specifications	Content	Application
L-Ornithine L-Aspartate	National Institutes for Food and Drug Control	50 mg/bottle	96.6%	For identification and content determination



3.2.S.7 안정성

3.2.S.7.1 안정성 요약과 결론

엘-오르니틴-엘-아스파르트산은 장기보존 안정성시험 및 가속 안정성시험을 실시하였다.

1) 실시한 안정성 시험 조건

Division	Condition	Period	Methods
Long-term stability	25±2℃, 60±10%RH	0,3,6,9,12, 18, 24,36 month	DAB
Accelerated stability	40±2℃, 75±5%RH	0,1,2,3,6 month	DAB

*According to ICH guideline ICH Topic Q1A(R2), ICH Topic Q1B)

Test items

content	Test items
Long-term stability	Appearance、 Water、 Clarity / Coloration、 Related substance、 Assay
Accelerated stability	Appearance、 Water、 Clarity / Coloration、 Related substance、 Assay

2) 사용한 Batch 정보

Batch Number	Batch Size	Manufacture Date	Expiration Date
C552202010	400.84kg	2022.02.08	2024.02.07
C552202011	400.84kg	2022.02.09	2024.02.08
C552202012	400.84kg	2022.02.09	2024.02.08

4) 안정성시험 항목과 시험방법

- 안정성시험 방법은 3.2.S.4에 따라 진행하였다.

5) 안정성시험 실시 결과

PV 로 제조된 LOA 3 배치(C5522010, C552202011, C552202012)는 3 개월동안 가속안정성과 장기 안정성 테스트를 실시하였으며, 현재까지 유의미한 변화는 관찰되지 않았다.

한국의 경우 유효기간 3 년이 지났으며, 정상적인 저장 조건에서 엘-오르니틴-엘 -아스파르트산은 상온 $\leq 30^{\circ}\text{C}$ 에서 안정하다.

SUNGWUN PHARMACOPIA CO., LTD.

Test		Specification	Initial	1 M	3 M	6M
			2022.02.19	2022.03.21		
Appearance		White crystalline powder or colorless crystals and well soluble in water	Conform	Conform		
Identification	A. Specific rotation	+26.5 ~ +29.0 °	28.3	--		
	B. IR	Corresponding to STD	Conform	--		
Purity	1) Clarity/Coloration	Clear, Colourless	Conform	Conform		
	2) pH	6.0 ~ 7.0	6.3	6.2		
	3) Related substance (TLC)	Spots other than the main spots in the test solution are not larger than the spots obtained in the standard solution	Conform	Conform		
	4) Chloride	≤ 300 ppm	<300	--		
	5) Sulfate	≤ 200 ppm	< 200	--		
	6) Ammonium	≤ 400 ppm	<400	--		
	7) Iron	≤ 30 ppm	< 30	--		
	8) Heavy metals	≤ 10 ppm	< 10	--		
	9) Water	≤ 7.0 %	0.66	0.81		
	10) Sulfated ash	≤ 0.2 %	0.03	--		
Assay		98.0 ~ 102.0 %	98.9	98.6		
Residual solvent (Methanol, In-house)		≤ 3,000 ppm	95	--		

shelf of life of L-Ornithine-L-Aspartate: 3years from the manufacturing date.
Accelerated term test's all test are observed no significant changes.

1.3) C552202012: 40±2°C, 75±5%RH(Accelerated)

Test		Specification	Initial	1 M	3 M	6M
			2022.02.19	2022.03.21		
Appearance		White crystalline powder or colorless crystals and well soluble in water	Conform	Conform		
Identification	A. Specific rotation	+26.5 ~ +29.0 °	28.3	--		
	B. IR	Corresponding to STD	Conform	--		
Purity	1) Clarity/Coloration	Clear, Colourless	Conform	Conform		
	2) pH	6.0 ~ 7.0	6.3	6.3		
	3) Related substance (TLC)	Spots other than the main spots in the test solution are not larger than the spots obtained in the standard solution	Conform	Conform		
	4) Chloride	≤ 300 ppm	< 300	--		
	5) Sulfate	≤ 200 ppm	< 200	--		
	6) Ammonium	≤ 400 ppm	< 400	--		
	7) Iron	≤ 30 ppm	< 30	--		
	8) Heavy metals	≤ 10 ppm	< 10	--		
	9) Water	≤ 7.0 %	0.77	1.27		
	10) Sulfated ash	≤ 0.2 %	0.01	--		
Assay		98.0 ~ 102.0 %	99.0	98.6		
Residual solvent (Methanol, In-house)		≤ 3,000 ppm	89	--		
shelf of life of L-Ornithine-L-Aspartate: 3years from the manufacturing date. Accelerated term test's all test are observed no significant changes.						

Sung

Long-term test shall be conducted for 36 months in accordance with the stability test plan.

2.3) C552202012: 30±2℃, 65±5%RH(long-term)

Test		Specification	Initial 2022.02.11	3 M 2022.05.20	6 M	9 M	12 M	18M	24 M	36 M
Appearance		White crystalline powder or colorless crystals and well soluble in water	Conform	Conform						
Identification	A. Specific rotation	+26.5 ~ +29.0 °	28.3	--						
	B. IR	Corresponding to STD	Conform	--						
Purity	1) Clarity/Coloration	Clear, Colourless	Conform	Conform						
	2) pH	6.0 ~ 7.0	6.3	--						
	3) Related substance (TLC)	Spots other than the main spots in the test solution are not larger than the spots obtained in the standard solution	Conform	Conform						
	4) Chloride	≤ 300 ppm	<300	--						
	5) Sulfate	≤ 200 ppm	< 200	--						
	6) Ammonium	≤ 400 ppm	<400	--						
	7) Iron	≤ 30 ppm	< 30	--						
	8) Heavy metals	≤ 10 ppm	< 10	--						
	9) Water	≤ 7.0 %	0.77	1.16						
	10) Sulfated ash	≤ 0.2 %	0.01	--						
Assay		98.0 ~ 102.0 %	99.0	99.2						
Residual solvent (Methanol, In-house)		≤ 3,000 ppm	89	--						

shelf of life of L-Ornithine-L-Aspartate: 3years from the manufacturing date.
 Long-term test shall be conducted for 36 months in accordance with the stability test plan.